



## **Terminology:**

- **Minimum level (ML)** represent the lowest quantifiable concentration in a sample based on the proper application of method-specific analytical procedures and the absence of matrix interferences.
- **Reporting limits (RL)** used by a laboratory may be from the analytical method, derived specifically for that laboratory, instrument or sample matrix; or as directed by a regulatory agency.
- **Method Detection Limit (MDL)** based on the ability of a measurement method to detect an analyte in the absence of a matrix (although sometimes in a sample matrix)

not detected (ND)	detected, not quantifiable (DNQ)	acceptable values		
	Method Detection Limit	Reporting Limit (RL) Minimum Level		



## Acceptable values

If your result is greater than or equal to the ML or RL:

- a) Use '=' as the Qualifier
- b) Enter the test result
- c) Ask your Regional Board contact if the MDL, ML or RL is required

	F	G	Н	1	J	K	L	Μ	N	ſ
1	Collection Date	Collection Time	Analysis Date	Qualifier	Result	Units	MDL	ML	RL	
125	8/15/2011	2:00	8/15/2011	=	1	mg/L				
126	8/16/2011	3:00	8/16/2011	=	• 0.6	mg/L				
127	8/17/2011	=			1.6	mg/L				
128	8/18/2011	<			1.4	mg/L				
129	8/19/2011	>			1.2	mg/L				
130	8/20/2011	>=			1.3	mg/L				
131	8/3/2011	ND			1.1	ug/L				
132	8/4/2011	12.30	0/4/2011	_	1.2	ug/L				



## Detected, not quantifiable (DNQ)

If your result is less than ML or RL but greater than or equal to the MDL:

- a) Use 'DNQ' as the Qualifier
- b) Enter the test result
- c) Enter the **MDL** and [**ML or RL**] along with a QA code if available

	F	G	Н		J	K	L	Μ	N
1	Collection Date	Collection Time	Analysis Date	Qualifier	Result	Units	MDL	ML	RL
125	8/15/2011	2:00	8/15/2011	=	1	mg/L			
126	8/16/2011	3:00	8/16/2011	DNQ	• 0.3	mg/L	0.1	0.5	
127	8/17/2011	=			1.6	mg/L			
128	8/18/2011	<			1.4	mg/L			
129	8/19/2011	>			1.2	mg/L			
130	8/20/2011	>=			1.3	mg/L			
131	8/3/2011	ND			1.1	ug/L			
132	8/4/2011 <sup>I</sup>	12.30	0/4/2011	_	1.2	ug/L			
		· · -			· -	4.			



## Not detected (ND)

If the analytical result is less than the MDL:

- a) Use 'ND' as the Qualifier
- b) Leave the result blank
- c) Enter the MDL

	F	G	Н		J	K	L	М	Ν
1	Collection Date	Collection Time	Analysis Date	Qualifier	Result	Units	MDL	ML	RL
125	8/15/2011	2:00	8/15/2011	=	1	mg/L			
126	8/16/2011	3:00	8/16/2011	ND	<b>•</b>	mg/L	0.1		
127	8/17/2011	=			1.6	mg/L			
128	8/18/2011	<			1.4	mg/L			
129	8/19/2011	>			1.2	mg/L			
130	8/20/2011	>=			1.3	mg/L			
131	8/3/2011	ND			1.1	ug/L			
132	8/4/2011 <sup>l</sup>	12.30	0/4/2011	-	1.2	ug/L			